Das Herbst®-Scharnier

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CE 0483

DENTAURUM

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The Herbst Appliance, Hans Pancherz, Editorial Aguiram

ISBN: 84-88453-07-4

Pancherz, H.: Treatment of class II malocclusions by jumping the bite with the Herbst-appliance; Am. J. of Orth. **76**, 343-348, 1980.

Weisländer, L.: Long-term effect of treatment with the headgear-Herbst appliance in the early mixed dentition. Stability or relapse?, Am. J. Orthod. Oct; **104**: 319-29 1993.

Kucukkeles, N. and Sandalli, T: Cephalometric evaluation of the therapeutic effects of the Herbst appliance in the treatment of class II. Div. I malocclusion, J. Marmara Univ. Dent Fac Sept; **1**: 230-6, 1992.

Windmiller, EC: the acrylic – splint Herbst appliance: a cephalometric evaluation., Am J Orthod. Jul; **104**: 73-84 1993.

Pancherz, H. and Anehus-Pancherz, M: The headgear effect of the Herbst appliance: a cephalometric long-term study., Am J Orthod. Dentofacial Orthop. Jun; **103**: 510-20 1993.

Pangrazio-Kulbersh, V and Berger, J. L.: Treatment of identical twins with Frankel and Herbst appliances: a comparison of results, Am J Orthod. Dentofacial Orthop. Feb; **103**: 131-7 1993.

Hansen, K: Post-treatment effects of the Herbst appliance. A radiographic clinical and biometric investigation, Swed Dent J Suppl; **88**: 1-49 1992.

Hagg, U: Change in mandibular growth direction by means of a Herbst appliance? A case report, Am J Orthod. Dentofacial Orthop. Nov; **102**: 456-63 1992.

Schiavoni, R. Grenga, V and Macri, V: Treatment of Class II high angle malocclusions with the Herbst appliance: a cephalometric investigation., Am J Orthod. Dentofacial Orthop. Nov; **102**: 393-409 1992.

Hansen, K and Pancherz, H: Long term effects of Herbst treatment in relation to normal growth development: a cephalometric study., Eur J. Orthod. Aug; **14**: 285-95 1995.

Hansen, K and Pancherz, H and Haag, V: Long-term effects of the Herbst appliance in relation to the treatment growth period: a cephalometric study., Eur J. Orthod. Dec; **13**: 471-84 1991

Pancherz, H: The Nature of Class II relapse after Herbst appliance treatment: a cephalometric long-term investigation., Am J Orthod. Dentofacial Orthop. Sep; **100**: 220-33 1991.

Jack, A: The Herbst appliance., Dent Tech. Mar; **44**: 10-1,15-6 1991.

Hansen, K, Pancherz, H and Peterson, A: Long-term effects of the Herbst appliance on the cranio mandibular system with the special reference to the TMJ., Eur J Orthod. Aug; **12**: 244-53 1990.

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- McNamara, J. A., Jr., Howe, R. P. and Dischinger, T. G.: A comparison of the Herbst and Frankel appliances in the treatment of Class II malocclusion., *Am J Orthod. Dentofacial Orthop.* Aug; **98**: 134-44 1990.
- Pancherz, H. and Fackel, O.: The skeletofacial growth pattern pre- and post-dentofacial orthopaedics. A long-term study of Class II malocclusions treated with the Herbst appliance., *Eur J Orthod.* May; **12**: 209-18 1990.
- Valant, J. R.: increasing maxillary arch length with a modified Herbst appliance., *J Clin Orthod.* Dec; **23**: 810-4 1989.
- Champagne, M: Herbst appliance therapy related to the mandibular plane angle., *Funct Orthod.* Nov-Dec; **6**: 17-21 1989.
- Dischinger, T. G.: Edgewise bioprogressive Herbst appliance., *J Clinical Orthod.* Sep; **23**: 608-17 1989.
- Kaplan, R. G: included condylar growth in a patient with hemifacial microsomia., *angle ortho.*, Summer; **58**: 85-90 1989.
- Bakke, M and Paulsen, H. U.: Herbst treatment in late adolescence: clinical, electromyographic, kinesiographic, and radiographic analysis of one case. *Eur J. Orthod.* Nov; **11**: 397-407 1989.
- Pancherz, H, Malmgren, O., Hagg, U., Omblus J. and Hansen, K: Class II correction in Herbst and Bass therapy , *Eur J orthod.* Feb; **11**: 17-30 1989.
- Valant, J. R. and Sinclair, P. M.: Treatment effects of the Herbst appliance, *Am J Orthod. Dentofacial Orthop.* Feb; **95**: 138-49 1988.
- Owen, A. H, 3d: Unexpected TMJ responses to functional jaw orthopedic therapy, *Am J orthod Dentofacial Orthop.* Oct; **94**: 338-49 1988.
- McNamara, J. A. and Howe, R. P.: Clinical management of the acrylic splint Herbst appliance, *Am J Orthod. Dentofacial Orthop.* Aug; **94**: 142-9 1988.
- Rider, E. A: Removable Herbst appliance for treatment of obstructive sleep apnea.; *J Clin Orthod.* Apr; **22**: 256-7 1988.
- Williamson, E. H: Skeletal change with long term wear of the removable Herbst in a post-adolescent male with anterior disk displacement., *Facial Orthop. Temporomandibular Arthrol.* Feb; **5**: 3-5 1988.
- Schiavoni, R: A method used to anchor the Herbst appliance.; *Cranio.* Jul; **6**: 245-51 1981.
- Hagg, U and Pancherz, H: Dentofacial orthopaedics in relation to chronological age, growth period and skeletal development. An analysis of 72 male patients with Class II division I malocclusion treated with the Herbst appliance., *Eur J. Orthod.* Aug; **10**: 169-76 1988.
- Pancherz, H and Hansen, K: Mandibular anchorage in Herbst treatment, *Eur J Orthod.* May; **10**: 149-64 1988.
- McNamara, J. A: Fabrication of the acrylic splint, Herbst appliance., *Am J Orthod Dento-facial Orthop.* Jul; **94**: 10-8 1988.

Handelman, C. S: Lower premolar extraction combined with Herbst treatment for mandibular retrognathia: a risky concept (letter), *Am J orthod Dentofacial orthop.* May; **93**: 445-6 1988.

Riolo, M. L., Moyers, R. E. and Tenhave, J. R: Imprecision and bias in orthodontic treatment results., *Am J Orthod. Dentofacial orthop.* Feb; **93**: 138-42 1988.

Howe, R. P: Lower premolar extraction/ removable plastic Herbst treatment for mandibular retrognathia., *Am J Orthod Dentofacial Orthop.* Oct; **92**: 275-85 1987.

Howe, R. P: Removable plastic Herbst retainer., *J Clin Orthod.* Aug; **21**: 533-7 1987.

Woodside, D. G., Metaxas, A. and Altuna, G.: The influence of functional appliance therapy on glenoid fossa remodelling., *Am J Orthod. Dentofacial Orthop.* Sep; **92**: 181-98 1987.

Pancherz, H and Hansen, K.: Occlusal changes during and after Herbst treatment: a cephalometric investigation, *Eur J Orthod.* Nov; **8**: 215-28 1986.

Williamson, E. H.: Treatment of anterior disk displacement with the removable Herbst appliance, *Facial Orthop Temporomandibular Arthrol.* Dec; **2**: 4-6 1985.

Goodman, P. and McKenna, P: Modified Herbst appliance for the mixed dentition.; *J Clin orthod.* Nov; **19**: 811-4 1985.

Williamson, E. H.: Removable Herbst treatment of anterior disk dislocation, *Facial Orthop Temporomandibular Arthrol.* Nov; **2**: 12-5 1985.

Pancherz, H. and Hagg, U: Dentofacial orthopedics in relation to somatic maturation. An analysis of 70 consecutive cases treated with the Herbst appliance, *Am J Orthod.* Oct; **88**: 273-87 1985.

Williamson, E. H.: Removable Herbst appliance utilization in adult disk displacement., *Facial Orthop Temporomandibular Arthrol.*, Jun; **2**: 11-2 1985.

Pancherz, H: The Herbst appliance – its biologic effects and clinical use., *Am J Orthod.* Jan; **87**: 1-20 1985.

Breads, P. R and Abruzzo, M: The Herbst appliance., *Quintessence Dent Technol.* Apr; **8**: 249-53 1984.

Howe, R. P.: The acrylic-splint Herbst, Problem solving., *J Clin Orthod.* Jul; **18**: 497-501 1984.

Howe, R. P. and McNamara, J. A. Jr.: Clinical management of the bonded Herbst appliance, *J Clin Orthod.* Jul; **17**: 456-63 1983.

Howe, R. P. Updating the bonded Herbst appliance., *J Clinical Orthod.* Feb; **17**: 122-4 1983.

Wieslander, L.: JCO interviews Dr. Lennart Wieslander on Dentofacial orthopedics. Head-gear – Herbst treatment in the mixed dentition, *J Clin orthod.* Aug; **18**: 551-64 1984.

Pancherz, H.: Vertical dentofacial changes during appliance treatment. A cephalometric investigation.; *Swed Dent J. Suppl.* **15**: 189-96 1.

Pancherz, H.: The mechanism of Class II correction in Herbst appliance treatment. A cephalometric investigation., *Am J orthod.* Aug; **82**: 104-13 1982.

Sarnas, H.: Pancherz, H., Rune B and Selvik, G.: Hemifacial microsomia treated with the Herbst appliance. Report of a case analyzed by means of roentgen stereometry and metallic implants, Am J orthod. Jul; **82**: 68-74 1982.

Howe, R. P: The bonded Herbst appliance, J Clin Orthod. Oct, **16**: 663-7 1982.

Langford, N. M, Jr.: Updating fabrication of the Herbst appliance, J Clin Orthod. Mar; **16**: 173-4 1982.

Pancherz, H and Anehus-Pancherz, M.: The effect of continuous bite jumping with the Herbst appliance on the masticatory system: a functional analysis of treated class II malocclusions, Eur J Orthod. Feb; **4**: 37-44 1982.

Langford, N. M, Jr.: The Herbst appliance, J Clin Orthod. Aug; **15**: 558-64 1981.

Pancherz, H.: The effect of continuous bite jumping on the dentofacial complex: a follow-up study after Herbst appliance treatment of class II malocclusions, Eur J Orthod. **3**: 549-60 1981.

Pancherz, H. and Anehus-Pancherz, M.: Muscle activity in Class II, division I malocclusions treated by bite jumping with the Herbst appliance. An electromyographic study, Am J Orthod. Sep; **78**: 321-9 1980.

Pancherz, H.: Treatment of Class II malocclusions by jumping the bite with Herbst appliance. A cephalometric investigation; Am J Orthod. Oct; **76**: 423-42 1979.

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